

# INFUSION PUMP STAND | MK-IS06

Zhangjiagang Medi Medical Equipment Supplies solution for hospital entire projects.  
including the hospital furniture and medical equipments



## PRODUCT DETAILS

### PUSH HAND



#### PUSH HAND

ABS push hand for easy movement

### HOOK



#### HOOK

Stainless steel four hooks, each hook can bear 2.5kg

### MOUNT BRACKET



#### MOUNT BRACKET

Mount bracket board by aluminum alloy material, accept customize drilling

### BASE



#### BASE

Aluminum alloy base by casting technique production, stable and firm

### CASTORS



#### CASTORS

3-inch silence medical casters, three with brake

### BASKET



#### BASKET

Basket for handle monitor/ECG wire accessories



## MAIN FEATURES

- Over all size: 610\*580\*1580-1660mm
- Height adjustment range 1580-1660mm by manual knob
- This cart can combine infusion pump and IV drip to use
- High grade stainless steel lifting column for fixing medical pumps
- Loading weight 60 lbs for column lift
- ABS push hand for easy movement
- Basket for handle monitor/ECG wire accessories

BETTER  
CARE  
BETTER  
COMFORT



**medik**

Tel : +86 15250358182

Fax : +86 512 58688450

E-mail : [info@medik-medical.com](mailto:info@medik-medical.com)

Address : Do.3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng  
Town, Zhangjiagang, Jiangsu Province, China



SCAN CODE