

PATIENT MONITOR STAND | MK-IS07

Zhangjiagang Medi Medical Equipment Supplies solution for hospital entire projects.
including the hospital furniture and medical equipments



PRODUCT DETAILS

PUSH HAND



PUSH HAND

ABS push hand for easy movement

MOUNT BRACKET



MOUNT BRACKET

Mount bracket board by aluminum alloy material, accept customize drilling

HEIGHT ADJUSTMENT



HEIGHT ADJUSTMENT

Height adjustment range 1340-1820mm by manual switch

BASE



BASE

Aluminum alloy base by casting technique production, stable and firm

CASTORS



CASTORS

3-inch silence medical casters, three with brake

BASKET



BASKET

Basket for handle monitor/ECG wire accessories



MAIN FEATURES

- Over all size: 650*650*920-1110mm
- Height adjustment range 920-1110mm by manual knob
- Aluminum alloy base by casting technique production, base size 650mm
- High grade stainless steel lifting column for fixing medical pumps
- 3-inch silence medical casters, three with locking
- ABS push hand for easy movement
- Basket for handle monitor/ECG wire accessories

BETTER
CARE
BETTER
COMFORT



medik

Tel : +86 15250358182

Fax : +86 512 58688450

E-mail : info@medik-medical.com

Address : Do.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng
Town,Zhangjiagang,Jiangsu Province,China



SCAN CODE