




## MAIN FEATURES

- \* Bed Frame Electrophoresis and powder coating
- \* Bed platform MS perforated platform with ventilation hole, molded
- \* Motor Quiet and robust electric actuators provide reliable operation
- \* Hand remote controller of all electrically operated functions and movements
- \* Handrails Collapsible aluminum railing with 6 bars
- \* Head and Foot board ABS material, detachable
- \* Castors: Caster diameter of 125mm with diagonal locking
- \* Bumper castors Bumper castors on 4 corners of bed
- \* Auto-counter function Auto-counter to provide the patient from sliding onto foot end


**medik**  
 Tel : +86 15250358182  
 E-mail : info@medik-medical.com




**YA-D3-3**  
 ELECTRIC HOSPITAL BED



**BACK SECTION**



**FOOT SECTION**



**HEIGHT ADJUSTMENT**

## PRODUCT DETAILS



## TECHNICAL PARAMETER

ITEM	PARAMETER
Overall size	2130*970*550-750mm
Lying size	1880*830mm
Back section up/down	0-75°
Foot section up/down	0-40°
High-low	550-750mm
Safe loading capacity	180kg
Lying board material	ABS plastic
Autoregress function	100mm
Castor	Dia125mm
Motor Control System	Jiechang or Linak or others