



## MAIN FEATURES

- \* Four split side rail, drop down
- \* Incorporated railing control with easy identifiable icon
- \* Nurse control incorporated in foot panel with LED display for moving state
- \* Detachable head/foot board with easy lock
- \* Central brake system with wheel diameter of 125mm
- \* Auto-regress to remove the pressure of the abdomen when adjust the backrest
- \* Auto-counter to provide the patient from sliding onto foot end
- \* Full length X-ray translucent platform
- \* Weighing scale system
- \* One button emergency stop switch can quickly stop all going movements.
- \* Electric CPR
- \* Manual CPR on both sides
- \* Bumper caster on four corners
- \* Angle Indicator for Backrest & Trendelenburg & Rev.Trend.
- \* Electric power height adjustment with foot control (optional)

**medik**  
 Tel : +86 15250358182  
 E-mail : info@medik-medical.com



**YA-D7-1**  
 ELECTRIC ICU BED



## PRODUCT DETAILS



## TECHNICAL PARAMETER

ITEM	PARAMETER
Overall size	2280*1050*520-770mm
Back section adjustable	0-75°±5°
Foot section adjustable	0-40°±5°
Height adjustable	520-770mm
Trendelenburg	0-12°
Reverse-Trendelenburg	0-12°
Castor	Dia125mm
Loading capacity	230kg