

medik
WWW.MEDI-MEDICAL.COM

OVERBED TABLE

MODEL:YA-T01



CONVENIENT | SECURITY | STRONG | COMFORTABLE

B
E
T
T
E
R
C
A
R
E
.
B
E
T
T
E
R
C
O
M
F
O
R
T

PRODUCT DETAILS

BETTER CARE BETTER COMFORT



» Material

HPL table top
Aluminum profile for the frame, heavy duty and easy to operate

» Air spring for tilting

Air spring for tilting adjustable
360° rotation angle



» Book / Laptop rest

Book rest included with tilting table top
Size:665*420mm

» Air spring for Hi-low

Air spring for height adjustment
745-1000mm



» Caster

Low profile twin wheel castors allow for
improved underbed access

» TECHNICAL PARAMETER

ITEM	PARAMETER
Model	YA-T01
Overall size (L*W*H)	885*420*745-1000mm
Table top size	885*420mm
Tilting table top size	665*420mm
Thickness of table top board	10mm
Castor	4 x Dia65mm (2 w/brake)
Safe working load	45kg

» PRODUCT SHOW

MAIN FEATURES

- Phenolic compact laminate (HPL) table top
- Aluminum undercarriage facilitate easy cleaning
- Air control height adjustment lifting devices.
- Single hand variable height operation for ease of use
- Book rest included with tilting table top



USER-FRIENDLY DESIGN

One hand lifting for height adjustment
Without press the lever, specially designed to avoid the table falling incident when it is pushed up without consciousness



ZHANGJIAGANG
MEDI MEDICAL EQUIPMENT CO.,LTD



Tel: +86 15250358182

Fax: +86 512 5868 8450

E-mail: info@medik-medical.com

Address: NO. 3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng
Town, Zhangjiagang, Jiangsu Province, China