

# INFUSION PUMP STAND | MK-IS05

Zhangjiagang Medi Medical Equipment Supplies solution for hospital entire projects.  
including the hospital furniture and medical equipments



## PRODUCT DETAILS

### MATERIAL



#### MATERIAL

Stainless steel body and aluminum alloy base, more durable

### HOOK



#### HOOK

Stainless steel four hooks, each hook can bear 2.5kg

### HEIGHT ADJUSTMENT



#### HEIGHT ADJUSTMENT

Height adjustment range 1340-1820mm by manual switch

### BASE



#### BASE

Aluminum alloy base by casting technique production, stable and firm

### CASTORS



#### CASTORS

3-inch silence medical casters, three with brake

### OPTIONAL BASKET



#### OPTIONAL BASKET

Basket for handle monitor/ECG wire accessories



## MAIN FEATURES

- Over all size: 650\*650\*1340-1820mm
- Height adjustment range 1340-1820mm by manual switch
- Aluminum alloy base by casting technique production, base size 650mm
- High grade stainless steel lifting column for fixing medical pumps
- 3-inch silence medical casters, three with locking
- Loading weight 60 lbs for column lift
- Available with two or four detachable hooks

BETTER  
CARE  
BETTER  
COMFORT



**medik**

Tel : +86 15250358182

Fax : +86 512 58688450

E-mail : [info@medik-medical.com](mailto:info@medik-medical.com)

Address : Do.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng  
Town,Zhangjiagang,Jiangsu Province,China



SCAN CODE