


MAIN FEATURES

- * ABS head and foot board with holes for moving the bed easily
- * Bumper castors on 4 corners for protection
- * 2 Aluminum alloy side rails, collapsible(Fold down by press the red button)
- * Gas spring for back section adjustment and CPR function
- * Castors with central braking system and Cranks for movements
- * Telescopic I.V. pole with adjustable height
- * Drainage hooks at the side of bed



PRODUCT DETAILS



medik
Tel : +86 15250358182
E-mail : info@medik-medical.com



YA-H5-3
MANUAL HOSPITAL BED



TECHNICAL PARAMETER

ITEM	PARAMETER
Overall size	2160*960*500-750 mm
Back section up/down	0-75°
Foot section up/down	0-40°
High-low	500-750mm
Trendelenburg	0-12°
Reverse-Trendelenburg	0-12°
Safe loading capacity	175kg