





YA-DS-R01

RECLINING PHLEBOTOMY CHAIR



PRODUCT DETAILS



Head rest

Comfortable built-in headrest



Retractable lock bar

Retractable lock bar secures the chair in the desired position



Armrest

Comfortable, padded armrests with guards to keep upholstery free from damage.



Casters

5" Total lock casters allow for better handling and control.



Side panel

Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair



Side panel Locker

Side panels "pop off" for easy cleaning and maintenance



- Ideal for use in dialysis, blood collection, respiratory and cardiac care
- Three positions: upright, deep recline, and elevated leg rest/footrest
- Retractable lock bar secures the chair in the desired position
- Padded push-handle for easy maneuverability.
- Comfortable, padded armrests with guards to keep upholstery free from damage.
- Blow-molded side tray includes recess for cup
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Comfortable built-in headrest
- 5" Total lock casters (in rear) allow for better handling and control.

1	Model	YA-DS-R01
2	Overall size(H*W*D)	1200 x 635 x 1040mm
3	Seat size(H*W*D)	530 x 482x 482mm
4	Tray size	650*390mm
5	Floor to Seat Height	660mm
6	Width between arms	508mm
7	Loading capacity	113kg