

MC-D08

ELECTRIC GYNECOLOGY CHAIR



Medi Medical Equipment Co., Ltd.
Better care Better comfort

ELECTRIC GYNECOLOGY CHAIR

+ CE ISO13485

Dimension	Table top dimension: L 2000mm * W 600mm Chair size: L 1450mm * W 600mm Filth basin size: 5L
Structure	Epoxy coated steel frame, the table top is upholstered with PU foam
Loading capacity	180kgs
Technical data	1.HI-LOW:580-1000mm 0 2.Back section fold up $\geq 70^\circ$ and fold down $\geq 20^\circ$ 3.Seat section fold up angle(Trendelenburg) $\geq 40^\circ$ 4.The leg section can be removable The leg support holder can rotate and adjust in height
Accessories	Leg support holder: 1 pair Waist grab holder: 1 pair Paper roll holder 1 piece Filth basin 1 piece Clamp 1 pair Foot pedal control 1 piece Arm rest 1 pair (Optional) Anesthesia screen holder 1 piece (Optional)

