

**PRODUCT
DETAILS**



Touch screen

LED touch control is used for operation, weighing system, instruction manual, patient card.



The 5th wheel

5th wheel steering function helps the operator to easily move the bed to the required space.




THREE COLUMN

Column structure design is to provide better safety and easy without complicated connection



Pushing handle

Ergonomic pushing handle foldable located outside head panel is easy for transporting the bed



Control panel

Adjust the status of the hospital bed according to different needs



casters

Dia 125mm castors. Central brake system is accessible from both sides

- Multiple functions can be achieved with this 3 column bed
- The left and right tilt function can be realized by the column
- One button total chair position is a special feature of this bed
- X-ray translucent back section
- Fifth wheel for direction
- Electric CPR and manual CPR on both sides
- Angle Indicator for Backrest & Trendelenburg & Rev. Trend.
- Weighing scale system (option)

1	Overall size(H*W*D)	2020-2230*1050*570-820mm
2	Back section adjustable angle	0-75°±5°
3	Knee section adjustable angle	0-35°±5°
4	HI-LOW	550-950mm
5	Trendelenburg	0-14°±1°
6	Reverse-Trendelenburg	0-14°±1°
7	Tilt left / right	0-12°±1°
8	Loading capacity	230kg

