




YA-DS-D09

DIALYSIS CHAIR



**PRODUCT
DETAILS**

Pillow
Removable pillow can be easily attached to head section for comfort.

Hand control
Hand control with visible icons for easy operation. CPR button is incorporated for emergency case.

Casters
Elegant design silent durable caster with individual brake

Armrest
The cushioned armrest can move with the backboard when the back-section is being adjusted

Footrest
Foldable footrest can give foot support during a semi fowler position and prevent the patient from sliding downwards.

IV Pole holder
infusion pole can be mounted space saving directly beside the chair.

- Three positions: chair position, folwer position, trendelenburg position
- Retractable lock bar secures the chair in the desired position
- Padded push-handle for easy maneuverability.
- Extra-wide padded armrests maximize comfort
- Durable epoxy coated welded steel frame is easy to clean and maintain
- Articulated headrest quickly adjusts to maximize comfort

1	Model	YA-DS-D09
2	Overall size	920*850*560mm
3	Seat size	570*500mm
4	Back section adjustment	-10°~75°
5	Leg section adjustment	-80°~0°
6	Trendelenburg position	0-12°
7	Loading capacity	200 KG

