

YAD-167

ROTATING NURSING BED

MAIN FEATURES



- 4 section bed surface with wire mesh design
- Hand rails to further aid transfer onto and off the bed
- Remote hand control for movements
- Profiled bed rotates through 90° to become chair
- Chair tilts forwards to aid transfer between sitting and standing
- One button chair rotation position
- Safety feature will not allow bed to tilt unless fully rotated
- Auto-regression function
- Four noiseless swivel castor, two with brake
- Can be set up to rotate either to the left or right
- Comes complete with foam mattress
- Easy clean, low maintenance
- Option Choice: motor brand, battery back up, head/foot board color, castor size

medik

Tel : +86 15250358182
E-mail : info@medik-medical.com



Castor

Total 4 pcs quality noiseless swivel castor dia 125mm for easy movement, 2 castors with independent locking can provide excellent stability



Head board

environment friendly beech wooden material, smooth surface with beautiful design and washable



Hand control

Top leading control system with remote hand control, easy operation for patient and care giver.



Overall size(H*W*D)	2050*970*480-770mm
Height adjustment	480-770mm
Back section adjustment	0-70°
Leg section adjustment	0-40°
Trendeleburg	0-12°
Revere- Trendelenburg	0-12°
One button chair rotation position	Y
Safety loading capacity	180kg