

INFUSION CHAIR

MK-F02

The infusion chair with dining table is mainly used in infusion rooms or as a medical recliner. It is equipped with storage baskets and infusion stand. All the designs greatly improve the comfort of patients.

- The frame is made of powder painted steel
- With environmental protection PVC leather
- A puller on the right side of chair to control the backrest and legrest simultaneously (gas spring operated) to 3 different positions
- With a stainless steel I.V rope with 4 or 2 hooks



ITEM	PARAMETERS
Folded size	73*84*105cm(folding)
Unfolded size	73*130cm(unfolding)
Structure	Epoxy coated steel
Back section adjustable	95°-120°-150°
Leg section adjustable	5°-60°
Packing	0*60*35cm=0.19M3 G.W.:36kgs
Safety loading capacity	150kg



Handle adjustment

The chair can be converted into a recliner (by pulling the side lever) quickly, along with automatic conversion into a chair when sitting up (by pressing the foot section). It is convenient for the patient to rest during the treatment.

Storage shelf

The infusion cart is equipped with a storage basket, which can hold water cups, bags, clothes and other items with a large capacity.



Arm rest

The wooden hand rest with ergonomic design greatly improves the patient's comfort.

Infusion stand

2 or 4 stainless steel hooks, the maximum weight that each hook can bear is 5KG.



Bottom structure

The bottom of the product adopts a triangular structure, which is extremely stable when the infusion chair is unfolded and placed on the ground. There is a non-slip plastic pad on the bottom to prevent displacement when the patient sits down.