

YA-D5-9

Electric Orthopedic Traction Bed



*Better care
Better comfort!*

YA-D5-9

Electric Orthopedic Traction Bed

DESCRIPTION

Orthopedic Traction Bed is designed for diagnosis, orthopedic surgeries, correction, traumatic treatment for the patient with disorders of the bones, joints, muscles, ligaments, tendons.

Multiple function and the accessories perform expanding traction to different parts of the vertebra and the extremities without moving the human body.

FEATURES

- Detachable head/foot board with easy lock
- Manual CPR on both sides(optional)
- Central brake system on four corner
- Hand remote control
- Four split side rail,drop down under the bed platform
- Leg and hand traction frame

STANDARD SPECIFICATION

Model	YA-D5-9
Overall size	2270*1060*450-725mm
Backrest adjustable	0-80°
Footrest adjustable	0-40°
Height adjustable	450-725mm
Trendelenburg	0-12°
Reverse trendelenburg	0-12°
Loading capacity	205kg



REMOTE CONTROL

The integrated remote control with clear icons allows nurses to quickly and clearly understand the sections they need to control.



FRAME

The overall frame is made of durable aluminum alloy pipe, the hexagon tube is convenient for stability. The frame can be quickly adjusted via unique switch wrench without tooling.



LEG TRACTION

Located at the end of the bed, it can be moved left and right through the knob, which is convenient for leg or body traction.



HANDLE

With two ergonomic handles, different positions can be controlled by the knob for different height patient, which is convenient for patients to take advantage of when it is inconvenient to get up.



PP SAFETY RAILING

Tuck-away PP safety railing with damper for slow dropping, special patented design to avoid accidental falling from patient's on bed operation



CENTROL LOCK SYSTEM

The central lock control system is located in the four corners, so that the nurse can lock or unlock the bed in every position. With red and green, easy to distinguish.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 199 6269 8550

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US