

YA-D5-3

Electric ICU Hospital Bed



*Better care
Better comfort!*

YA-D5-3

Electric ICU Hospital Bed



DESCRIPTION

YA-D5-3 is a total updated design full electric hospital bed, It equips with Timotion brand motor control system and widely used in ICU room due to stable quality. It will provide an economic solution for some ICU room due to the embedded side railing control/patient control and the remote nurse hand control.

STANDARD SPECIFICATION

Overall size	L2150*W980*H450mm
Timation control system	✓
Height adjustment range	420-750mm
Back board	0-75°
Knee board	0-35°
Height range	420-750mm
Tren/Reverse Tren	0-13°/ 0-13°
Loading weight	250kg

MANUAL CPR

The manual CPR board locate in the two side of the back board. The back board can be released to flat shortly.



ELECTRIC CPR

Electric CPR button locate in remote hand control and side railing control.



WHEELS AND LOCKING

Dia 125mm double side mute wheels, swivel type with central locking.



BUMP CASTOR

There are 2 bump castors in each corner, it will protect bed well.



EMBEDDED NURSE CONTROL

Nurse can adjust the bed for different position via the embedded button in the side railing.



PATIENT CONTROL

Patient can adjust the back and knee board only.



DEGREE INDICATOR

Nurse can adjust the bed more accuracy via the degree indicator.



DETACHABLE HEAD/ FOOT BOARD

Head/foot board can be removable easily and shortly due to the detachable design.



REMOTE HAND CONTROL

Nurse can adjust the bed via the remote hand control for different position.



DETACHABLE ABS PLATFORM

4 section detachable ABS platform, easily for clean.



POWERFUL BEAM SUPPORT

It will provide heavy loading support due to the special structure.



ELECTRIC CONTROL SYSTEM

Timotion control system.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 199 6269 8550

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US