

# MC-C14

Delivery Bed



*Better care  
Better comfort!*

# MC-C14

# DELIVERY BED

## DESCRIPTION

MC-C14 birthing bed provides a solution for pre-delivery- delivery- after-delivery on one bed, it can be converted from a general hospital bed to a delivery table without transferring the mother. The compact design with practical functionality makes the bed popular with limited budget hospitals while requiring LDR concept for birthing progress.

## FEATURES

- Three-section stop. Backrest is adjustable by gas spring.
- The leg section can be telescoped under the backrest for lithotomy positions.
- Removable Stainless steel head & foot boards with laminated panels.
- Foot support for high/low chair positions
- Complete with SS side railings, basin and telescopic IV pole
- Height is adjustable by crank.
- Four casters with diameter of 125mm, 2pcs with brakes

## STANDARD SPECIFICATION

Overall size	2000*750*500-750mm
Height adjustment	500-750m
Back section adjustment	0-70°
Loading capacity	220kg
Leg section adjustment	Pack up



**Back section**

Back section adjust by gas spring.



**Height adjust**

The height of delivery bed is adjusted by the crank, and the stroke is 250mm.



**Leg section**

Can be retracted into the bed to gynecological and delivery surgical.



**Mattress**

Flat and thick sponge mattress, more comfortable.



**Side rail**

Stainless steel side rail can be folded up and down.



**Castor**

Four 125mm diameter independent casters, 2pcs with brakes.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US