

medik

BETTER CARE · BETTER COMFORT

ANESTHESIA TROLLEY

MODEL:MK-C06



CONVENIENT | SECURITY | STRONG | COMFORTABLE

PRODUCT DETAILS

medik



» Tilt Storage Bin

Total size 620*340*110mm, five pcs anesthesia boxes each line, two lines

» Extended Work Surface

Easy for care giver writing



» Structure

SS 304# plate cover on top
Aluminium-plastic side panels for light weight cart with strength

» Sharp Container / Stainless Steel Basket



» Waste Bin

Two waste bins with lid

» 5" Braking Caster

Noiseless casters with heavy duty bearing design, 2 with brake



» Drawer

Electrolytic panels drawer for long-life usage

» TECHNICAL PARAMETER

ITEM	PARAMETER
Model	MK-C06
Overall size (L*W*H)	640*1900*420mm
Material	Steel/Aluminum/Polymer
Color	White,Red,Blue
Castor	Dia 125mm
G.W.	57.8KG + 20.4KG
Package size	680*530*1150mm + 880*200*670mm

» PRODUCT SHOW

MK-C06 Anesthesia Cart is a safe way of transporting medications, controlled substances and other supplies, from one department to another. Each cart comes with ample storage for needed materials with dividers and compartments available for improved organization. Medical Anesthesia Carts are ideal for all hospitals and surgery centers.





ZHANGJIAGANG
MEDI MEDICAL EQUIPMENT CO.,LTD



Tel: +86 15250358182

Fax: +86 512 5868 8450

E-mail: info@medik-medical.com

Address: NO. 3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng
Town, Zhangjiagang, Jiangsu Province, China