

MK-PC03

Medical Laptop Cart



*Better care
Better comfort!*

MK-PC03

MEDICA LAPTOP CART

DESCRIPTION

MK-PC03 is a new creative medical cart to utilize the platform with patient information management for dairy room check-up and practical medical stuff storage. With computer built in and various storage space in the cart, medical staff can improve the efficiency of workflow and optimize the hospital spaces.

FEATURES

- Height adjustable supporting column can be adaptable for medical staff in different heights
- Large working top table with integrated handle to move the cart easily
- Telescopic keypad and mouse tray
- Drawers and baskets for storage
- Mobile 4"(φ100mm) casters

STANDARD SPECIFICATION

Overall size	480 *520*938~1138mm
Working surface	480 *520mm
Height adjustable	938~1138mm



Slide-out Keyboard with Mouse Tray

The Slide-out Keyboard is easy to operate and has a place for the mouse below. The mouse tray can be on the left side or on the right side, which is convenient to meet different needs.



Upper Drawer

There is a drawer at the upper where the divider can be placed, which can improve utilization of the drawer.



Lower drawers

There are two drawers at the lower end to store sundries, and there are storage baskets on both sides, which greatly increases the storage space.



Height adjust

The height is controlled by a gas spring, and it can stay at the desired height when the pedal is released, meeting the needs of different heights.



Laptop lock (optional)

There is a Laptop lock on the desktop, with a key that can adjust the length to adapt to different sizes laptops, and fix the laptops to prevent slipping and theft.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US