



MODEL : YA-AC01
AMBULANCE STRETCHER

YA-AC01

Ambulance stretcher

The YA-AC01 is a tough, all-aluminium manual ambulance cot. versatile and stable, it has a retractable head section for 360-degree mobility in any height position. YA-AS01 is widely used in various emergency scenarios.

MAIN FEATURES

- One person can operate and send the patient to the ambulance vehicle
- The surface is made of foamed cushion, vinyl protective cover that is durable and easy to clean
- Equip I.V. pole
- Castors dia 150mm

TECHNICAL SPECIFICATION

Load capacity	180KG
Size(unexpanded)	206*57*28CM
Size(Expanded)	191*57*91CM
Material	High-strength aluminum alloy





Fold-down side rails

The foldable guardrail can prevent the patient from falling off the stretcher and causing secondary injury.



Legs folded crank

When pushing the stretcher to the ambulance, you can press the handle of the head to help retract the legs, which is quick and convenient to transfer the patient to the ambulance.



Large wheels

Large wheels with a diameter of 150mm, with non-slip patterns, can easily cross various obstacles, and with a lock system, which can fix the stretcher in the specified position.



Back adjust

The back is controlled by a gas spring, and the back can be adjusted to the desired angle by just pressing lightly.



Bolster mattress

The concave mattress provides support to the patient while minimizing the shaking of the patient and avoiding secondary injuries.



Auxiliary wheel

There are spare wheels on the four corners of the stretcher, which can be moved by the wheels even after the stretcher is pushed onto the ambulance, which saves effort.



You Tube



Medik

Tel : +86 13962231292 / Fax : +86 51258688450

E-mail : info@medik-medical.com

Address : No.3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng Town, Zhangjiagang, Jiangsu Province, China