



YA-AS07

DETACHABLE
STRETCHER

**BETTER CARE
BETTER COMFORT**

Zhangjiagang
Medi medical Equipment
Co., Ltd



Tel : +86 15250358182 / Fax : +86 51258688450
E-mail : info@medik-medical.com
Address : No.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng
Town,Zhangjiagang,Jiangsu Province,China



Height adjust

The stretcher has the function of up and down, and can be placed in any position to meet the separation of the bed surface, which is faster and more convenient.



Fixed slot

The stretcher has the function of up and down, and can be placed in any position to meet the separation of the bed surface, which is faster and more convenient.



IV hole

There are infusion holes on the four corners of the transfer vehicle to meet the needs of infusion.



Storage basket

There is a storage basket at the bottom of the stretcher to meet basic storage needs.



Back adjust

The back can present different angles to meet several needs.

A-AS07 is a detachable stretcher for ambulances

The bed surface can be quickly separated, suitable for transportation inside the hospital and emergency center transportation and outdoor rescue. With the function of height lifting and back adjustment, it can well meet the needs of various scenes.

Main features

- This multifunctional automatic loading stretcher has varied heights, it's made of high-strength aluminum alloy thick-wall pipes.
- Support for IV stand included
- Aluminum alloy frame, easy cleaning and sterilization.
- This type of stretcher can be characterized as light weight, durable, safe and easy to operate.

Technical data

ITEM	PARAMETERS
Overall size	190*55*54-88cm
Back adjust	0-75°
Load capacity	159KG

