

YA-PS10

Patient transfer stretcher



*Better care
Better comfort!*

YA-PS10

PATIENT TRANSFER STRETCHER

DESCRIPTION

The powder-coated structure, basic adjustment functions and smooth movement enable YA-PS10 to meet basic patient transfer and care needs, it is easy to clean and cost-effective. The storage space under the bed adds flexibility of use.

FEATURES

- Powder coating structure,easy to clean and durable
- Back section can be adjusted by gear
- Side rail can be raised and dropped down
- Mattress thickness 5cm
- Large storage basket
- Four castors with diagonal lock
- Knock down packing
- With IV stand

STANDARD SPECIFICATION

Overall Size	L1800*W500*H800mm
Side rail size	L700mm*H350mm
Mattress thickness	5cm
Packing	Knock down design
Safe Working Load	150kg
Castor	Dia 100mm



Back section adjustment

Back section adjusted by gear and can be stopped at any height.



Side rail

Side rail can be raised and dropped down.



IV stand

IV stand with removable hooks.



Storage space

Big storage basket at bottom of the stretcher.



Castors

4 units silent castors with diagonal lock, easy to move.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US