



YA-PS12

Medical Parallel Transport Bed



YA-PS12

MEDICAL PARALLEL TRANSPORT BED

Pain Points

High Transfer Difficulty	Risk Of Secondary Injury	Occupational Health Hazards	Medical Disputes
ICU patients, overweight patients, and post-operative patients present significant physical challenges for medical staff.	Improper manual transferring can cause bone dislocation, wound aggravation, or vascular tearing.	Medical workers frequently suffer from spinal and lumbar/sacral injuries due to heavy lifting.	Adverse medical consequences and patient discomfort during transfer often lead to legal and financial disputes.

Product Advantages



During transportation no relative movement between patient and the bed surface

The unique design of the device enables patients to safely, stably, and comfortably pass through different heights of surgical beds, hospital beds, and transfer beds in different positions, effectively alleviating the troubles of bed to bed transferring for those who with mobility difficulties.



Easy to operate

The product has an intelligent electric power assist system, adopts a high-quality lifting system, reliable braking, and flexible implementation. The intelligent control panel and electric drive move up and down, making it easy for family members or medical staff to operate alone. Complete safety precautions, convenient operation, time-saving and labor-saving, safe and fast.



Easy to clean and nursing for nurse

Nurse need to Clean patient and the hospital beds, this device great help to move the patient away shortly.



YA-PS12

MEDICAL PARALLEL TRANSPORT BED

Product Description

The gurney is used to care, move and transfer patients, the aged and the disabled, to complete automatic bed-crossing

1. To shifting patient to different position, allowing carer clean the patient and bed
2. To Transfer patient, including pick up patient from other bed or operating table, moving patient to new place.
3. One-man operation needed,with armature allow moving up and down, parallel moving, powered walking.
4. during moving /transferring, there is no relative movement between patient and the bed, comfortable and safety to patient.
5. Meet medical equipment basic standard. 6. With electric power steering and digital control is easy to save both time and labor, providing instant security.




YA-PS12

MEDICAL PARALLEL TRANSPORT BED

Functions



Integrated Console



Retractable Lying Board



Liftable Height



Movable Guardrail



Center Swivel Wheel

YA-PS12

MEDICAL PARALLEL TRANSPORT BED



Integrated Console Set

All-in-One Control Terminal

The ergonomic handle integrates a power switch, speed control, and direction toggles into a single terminal. A clear LCD screen provides real-time monitoring of speed and battery levels, ensuring a safe and effortless user experience.

Intelligent Power Assist

Designed to reduce the physical workload of caregivers, the electric assist system adapts to various slopes by adjusting torque output. Combined with a small turning radius, it allows for nimble and precise navigation in confined hospital spaces.



Center Swivel Wheel

Turn Freely

The retractable central swivel wheel offers dual modes via a simple pedal: retract it for stable, straight-line tracking, or lower it for 360° maneuverability, including turning in place.

Run Smoothly

Engineered for stability, the retractable wheel and large 200 mm casters absorb impact across uneven surfaces. This design ensures maximum patient comfort and effortless obstacle crossing.



Liftable Height

Versatile Application

The adjustable height allows seamless alignment with various hospital beds and diagnostic equipment. This flexibility ensures smooth patient transfers while significantly reducing the risk of spinal strain for caregivers.

Enhanced Safety

Lowering the vehicle to its minimum height provides a stable center of gravity during transit. This design facilitates emergency response and maximizes patient safety and comfort throughout the transfer process.

YA-PS12

MEDICAL PARALLEL TRANSPORT BED



Effortless Patient Transfer

Achieve seamless patient transfer between hospital beds and diagnostic equipment without manual lifting or third-party assistance. By utilizing the synchronized movement of the telescopic board and transfer cloth, this system eliminates the risk of secondary injuries and spinal strain, ensuring maximum safety for both patients and caregivers.



Movable Guardrail

Noise Reduction

The upgraded snap-on guardrails feature a one-click telescopic mechanism for smooth, quiet operation. This gentle process minimizes patient disturbance and eliminates the risk of hand-pinching for medical staff.

Smart Fall Prevention

Extended guardrails provide essential coverage of the patient's center of gravity to minimize fall risks. For maximum security, integrated sensors ensure the transfer system only activates when the guardrails are securely locked in position.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US