

# PATIENT TRANSFER STRETCHER

YA-PS14 is a manual patient transfer stretcher with height and back adjustment and exit-assisted guardrail. It can meet the general needs of hospital care. Equipped with 5th wheel can be moved easily.



## Back section

The back section is adjusted via gas spring to any wanted degree between 0~70



## Height adjustment by crank

The height adjustment is controlled by crank, which is convenient for nurses to control



## Aluminum alloy guardrail

The 4-piece aluminum alloy guardrail can maximize the length of the guardrail and prevent the patient from slipping off, meanwhile it can better assist the patient exit with half railing folded down.



## Central brake castor

Diameter 200mm big castor with central control brake has good passability during transportation



## 5th castor

With 5th retractable wheel in the middle for direction



## Oxygen cylinder holder

The oxygen cylinder holder is specially designed under the backboard of the stretcher frame



## IV Pole

The storage hole of the iv pole is hidden under the back board



- **Overall size**  
1930\*640\*520-820mm
- **Back section adjustment**  
0-70°
- **Height adjustment**  
520-820mm
- **Caster**  
Diameter 150mm
- **Loading Capacity**  
250kg

- Manual luxurious stretcher
- Side rails: 4pcs, Aluminum alloy material, foldable
- ABS platform easy for cleaning
- Manual height adjustable
- With 5th retractable wheel in the middle for direction
- Backrest adjustable in a pneumatic way
- Central brake system with 150mm diameter wheel



Zhangjiagang Medi Medical Equipment Co., Ltd  
Tel: +86 15250358182  
E-mail: info@medik-medical.com

