

YA-SS07

Lightweight Stair Stretcher



*Better care
Better comfort!*

YA-SS07

LIGHTWEIGHT STAIR STRETCHER

DESCRIPTION

The YA-SS07 is a light weight stair stretcher made of high strength aluminum alloy. It is widely used in the home, school, hospitals etc for those who can not move alone during emergency evacuation specially through stairs, which is very easy to store and transport with its light weight and folded structure.

FEATURES

- The stair stretcher is made of high-strength aluminum alloy materials.
- It is ideal for emergency staff to transfer patients up and down stairs.
- Under the front of the stretcher, there are two flexible uplift rods.
- The stretcher can be folded to save space in an ambulance or vehicle.
- It is lightweight, compact, easy to operate and also easy to sterilize and clean

STANDARD SPECIFICATION

| | |
|--------------------------|-------------------|
| Model | YA-SS07 |
| Size (Unfolded) | 900*510*910mm |
| Size (Folded) | 900*170*590mm |
| Loading capacity | 159kg |
| Carton Size (1pc/carton) | 960 * 160 * 540mm |
| N.W. | 8kg |
| G.W. | 9.2kg |



Material

Coated Oxford cloth is tough and not easy to damage, and the surface can be easily cleaned of blood stains and stains.



Castors

There are 4 castors which make it easier to transport the injured on a flat ground, and 2 with brakes, which can be fixed at designated positions.



Carrying handle

Adjustable carrying handles can be adjusted in length and angle for maximum convenience when going up and down stairs.



Folding handle

Stair stretcher can be folded to the smallest volume in one step for easy storage.



Foldable design

The whole stretcher include the handle can all be folded, which can reduce the volume and put it in any desired position.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun, Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.