

YA-CS01

Transfer Connecting Stretcher



*Better care
Better comfort!*

YA-CS01

TRANSFER CONNECTING STRETCHER

DESCRIPTION

YA-CS01 connecting stretcher is designed to minimize the contamination risk from transferring the outside stretcher into clean OT room directly. The upper stretcher can be sliding from the outside down trolley to the inside down trolley on wheels effortlessly, while each side stretchers is a normal stretcher for patient transfer. It is an ideal contamination control device for hospital.

FEATURES

- Two divided individual stretcher
- Mobile up stretcher sliding on aluminum orbit
- Central brake system
- Fifth wheel for direction
- Connecting system easy to operate
- Two long ABS side rail

STANDARD SPECIFICATION

Overall size	3650*640*640-970mm
Back section adjustment	0° to70°
Height adjustment	640-970mm
Castors	Four 150mm castors with a central locking system
Loading capacity	225kg



MAIN TROLLEY

The main trolley is made of ABS material, and the upper rails is made of aluminum alloy.



BACK SECTION

The back section is adjusted by gas spring to be between 0~75.



HEIGHT ADJUSTMENT

Height adjustment via crank from 640 to 970mm.



5TH DIRECTION WHEEL

The stretcher has 5th direction wheel provides the stretcher with ideal manipulation and stability.



CASTORS

The upper platform is transferred from one stretcher to another by means of casters.



CONNECT FOR STRETCHER

The two stretcher trolleys are connected by clips and the moving of the upper stretcher is realized.



CONNECTOR RELEASE

By pressing down the push handle of the main stretcher, the connecting stretcher can be separated.



IV POLE HOLDER

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US