

MC-C07

S.S. Gynecologist Chair With Drawer



*Better care
Better comfort!*

MC-C07

S.S. GYNECOLOGIST CHAIR WITH DRAWER

DESCRIPTION

The MC-C07 Gynecology Chair is specifically designed for use during a gynecological examination. This type of chair allows the patient to be comfortable and relaxed while the doctor performs the exam.

FEATURES

- The table is made in 3 sections.
- Framework made of Stainless Steel.
- With adjustable Backrest & leg rest to varying positions.
- Moulded seat with a waste tray underneath.
- Synthetic leather upholstery.
- Side step integrated on the side of the table.

TECHNICAL PARAMETER

Model	MC-C07
Overall size	0-75°
Back section adjustable	1870*600*680mm
Foot section adjustable	0-90°
Safe working load	150kg



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US