

MC-D02

Delivery Bed



*Better care
Better comfort!*

MC-D02

DELIVERY BED

DESCRIPTION

The MC-D02 bed is a comfortable maternity and delivery bed. The expecting, birthing and recovering woman can all make use of this bed. MC-D02 delivery bed is equipped with every essential feature and designed with all labor positions in mind.

FEATURES

- The base seat and height adjustment system are made of power coated steel, anti corrosion.
- The bed cover and head&foot board made by PP material, the head&foot board can be movable.
- Equipped with switch.
- The handle can be disassembled and moved.
- The leg support can be dismantled and movable.
- The Lithotomy rod can be dismantled.
- The mattress adopts high density memory sponge.
- Foot section can be retracted into the bed to gynecological and delivery surgical.

STANDARD SPECIFICATION

Overall size	2050*980*580-885mm
Back section adjustable	0-70°±5°
Height adjustable	580-885mm
Trendelenburg	12°
Reverse-trendelenburg	7°
Loading capacity	230kg



Bed Surface

The bed surface is composed of PP cover, easy to clean.




Guardrail

The guardrail can be rotated and lowered.



Mattress

The mattress adopts high density memory sponge and is once formed. Mattress can be removed.



Foot Section

Foot section can be retracted into the bed to gynecological and delivery surgical.



Emergency Stop Button

For emergency situations, all movements can be paused with one click.



Basin



Leg Support



Hand Controller



Castor



Handle

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US