

# MC-D09

Electric Birthing Bed



*Better care  
Better comfort!*

# MC-D09 ELECTRIC BIRTHING BED

## DESCRIPTION


The MC-D09 is a high-end birthing bed designed to support every stage of childbirth—from prenatal examinations to labor, delivery, and postpartum recovery. Engineered for comfort, safety, and clinical precision, the MC-D09 offers advanced adjustability, excellent mobility, and hygiene-focused materials to meet the demanding needs of modern obstetrics.

## FEATURES

- The bed design ensures safety of the mother while being on the bed and during labor.
- Three-section electric adjustment enables smooth control of height, backrest, and leg positions.
- Integrated design supports labor, delivery, and recovery in a single bed.
- Trendelenburg position allows flexible adjustments during labor or emergencies.
- Removable leg section allows quick transition between examination and delivery.
- Seamless antibacterial upholstery ensures comfort and easy cleaning.
- Handheld and panel controls provide precise, convenient operation.
- Central locking wheels ensure mobility and secure positioning.

Patient

- Freedom to change a birthing position
- Yes to natural childbirth
- Maximum safety for baby and woman in labor




- Intuitive and easy operation
- Hygiene without compromise
- Minimum physical effort



Staff

Management

- Exceptional modern design
- Increased level of care
- LDR and LDR Pconcept





**TECHNICAL PARAMETER**

Model	MC-D09
Total length of the lying area	2122 mm
Minimum length in the default setting	1570 mm
Full width	980 mm
Minimum height	600 mm
Maximum height	1050 mm
Back section up/down	70°/10°
Seat section up/down	18.5°/0°
Leg rests vertical	135°
Leg rests horizontal	60°
Foot section up/down	22°/ 0°
Trendelenburg position	10°
Castor size	150 mm
Maximum patient load	230 kg

**STANDARD ACCESSORIES**

Accessory latch	1 Pair
AISI 304 waste basin	1 Unit
Maternal assisting handle	1 Pair
Caster	4 Pc
Cable	1 Pc
Legrests	1 Pair
Mattress	1 Pc
IV pole socket	2 Pc
Foot switch (optional)	1 Unit
Backup battery (optional)	1 Set
S.S kicking bucket (Optional)	1 Pair
Hand remote controller (Optional)	1 Pair

**DETAILS****Midwifery handle**

Maternal assisting handle to provide assistance during childbirth.

**Manual CPR handle**

Manual CPR handle, to adjust backrest in flat position in emergency.

**AISI 304 accessory latch**

AISI 304 accessory latch, facilitate accessories installation.

**Castor**

Imported TENTE casters of diameter 150mm, with directional function.

**IV pole holder**

IV pole socket.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US