

MC-H04 [MANUAL OBSTETRIC TABLE]

MC-H04 type integrated obstetric bed is suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit.

Features

- ◆ The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.
- ◆ The main bed and auxiliary bed are connected with each other, auxiliary table is easy for dismantling.
- ◆ Foot step hydraulic elevating, Imported Y type sealing ring and is durable.
- ◆ Equipped with stainless steel sewage basin and other various accessories.
- ◆ The mattress adopts high density memory sponge and is once formed with no seam.
- ◆ Rotary brake device which is easy for moving obstetric table.



Standard Accessories: A pair of arm resting frame, a pair of leg support frame, a pair of armrest, a stainless steel sewage basin.

MC-H04A、MC-H04B [SPECIFICATION]

Mode	MC-H04A	Mode	MC-H04B
Table Length mm	1950±50	Table Length mm	1950±50
Table Width mm	600±10	Table Width mm	600±10
Height Adjustment mm	(750 ~ 950)±50	Height Adjustment mm	(750 ~ 950)±50
Leg board Upwards	50°	Backboard Upwards	50°
Leg board Downwards	10°	Backboard Downwards	10°
Seat board upwards	30°	Seat board upwards	30°
Seat board downwards	10°	Seat board downwards	10°
Working area	650 x 600	Working area	650 x 600

MC-H04C [SPECIFICATION]

Mode	MC-H04C
Table Length mm	1900±50
Table Width mm	600±10
Height Adjustment mm	(700 ~ 950)±50
Backboard Upwards	50°
Backboard Downwards	3°
Trendelenburg	10°
Rev-Trend	20°
Leg board Outwards	90°
Working area	600 x 600