



## MODEL : MC-C06B

SIMPLE GYNECOLOGICAL EXAMINATION CHAIR

## MC-C06B

MC-C06B is a Simple gynecological examination chair in knockdown design, using high-quality stainless steel plate and round tube with argon welding, smooth and bright, without concave deformation, can meet the needs of most gynecological examinations.

**1900\*610\*750mm**  
Overall size

**0-42°**  
Buttocks adjustment

**0-75°**  
Back section adjustment

**170kg**  
Loading capacity

**0-90°**  
Leg section adjustment



## Detailed description



### Back adjust

The back can be adjusted in multiple stages, using a ratchet structure.



### Buttocks adjustment

The buttocks can be lifted by Ratchet structure, which is convenient for the doctor to check.



### Leg support

The leg support can be adjusted in height and curvature, which greatly improves the comfort of the patient.



### Dirt Basin

Stainless steel dirt basin easy to clean, durable and replaceable.



### Bed legs

Rubber legs are installed on the bed legs to prevent noise from rubbing with the floor.



### Leg section adjustment

The leg section is also adjusted by the ratchet structure to 90 degrees down so the doctor can get closer with convenience, it also can be adjusted in a horizontal position with the seat section.



You Tube



Medik

Tel : +86 13962231292 / Fax : +86 51258688450  
E-mail : info@medik-medical.com  
Address : No.3 Jinxiu Road, Jinfeng Industrial Zone, Jinfeng Town, Zhangjiagang, Jiangsu Province, China