



MODEL: MC-D04B

GYNECOLOGICAL TABLE

GYNECOLOGICAL TABLE

MC-D04B gynecological table is mainly used for delivery and gynecological procedure. Equipped with hidden movable auxiliary units to facilitate the operation of medical staff, it is a good choice for hospital obstetrics and gynecology.

MAIN FEATURES

- All base and column cover is made of 304 medical grade stainless steel
- The ascent and descent, forward and backward tilting of this set are manipulated by micro touching controller, which is convenient and easy to operate.
- The back section adopts gas springs structure, and the angle can be regulated.
- Foot controls the fix and move of the base set.
- The basins and all outer shells are made of stainless steel, which makes nice appearance and easy cleaning.
- The choice of colors of foot and soft mat can be available.

TECHNICAL PARAMETER

ITEM	PARAMETER
Overall size	1840* 600 * 750-1000mm
Back section adjustment	0°-60°
Height adjustment	750-1000mm
Trendelenburg	0°-22°
Reverse trendelenburg	0°-22°
Loading capacity	200kg



ARMREST



HAND GRIP



HAND GRIP




HAND GRIP



ANESTHESIA SCREEN FRAME



* MC-D04B



FOOT SECTION

Can be retracted into the bed to gynecological and delivery surgical. Unfolding can become 3-section gynecological exam table



POWER SWITCH AND FOOT PEDAL

Control power, the motor will not work when turned off. Foot pedal control for height and trend and reverse trend, freeing medical staff's hands.



BACK SECTION

Back section is adjusted from 0° to 60° by gas spring

Zhangjiagang Medi Medical Equipment Co., Ltd.

Tel: +86 139 6223 1292
E-mail: info@medik-medical.com
Address: Room 1, No. 5, Xinzha East Road, Zhashang Village,
Yangshe Town, Zhangjiagang City, Suzhou, Jiangsu Province



You Tube

Medik