

MC-D06

Delivery Bed



*Better care
Better comfort!*

MC-D06 /Delivery Bed

MC-D06 is a fully electric LDR bed with closed structure and integrated guardrail control panel for easy clean and operation. The unique design of the leg section that fold up into footrests helps the patient adjust her delivery posture. The integrated design of the delivery bed can meet the needs of the patient during the entire delivery process needs, including labor, delivery and recovery stages, and improve physician efficiency and physician-patient comfort.

MAIN FEATURES

- Bed Frame made of steel with powder coating surface
- Closed structure, easy for clean
- With guardrail control panel
- Folded up design of the leg section, convenient for patient
- Removable stainless steel waste basin with drain pipe
- Three-position central brake castors with guidance function
- Trendelenburg and reverse trendelenburg position
- Two split side rails, drop down under the bed platform
- Detachable headboard with easy lock
- With angle indicator embedded in the railing
- Bumper castors at the head of the bed for protection



FUNCTIONS



Reverse trendelenburg

Reverse trendelenburg can be adjusted by remote control from 0-10°



Trendelenburg

Trendelenburg can be adjusted by remote control from 0~10°



Height adjustment

Height can be adjusted by remote control from 600~900mm



Back section adjustment

Back section can be adjusted by remote control from 0-65°

TECHNICAL PARAMETERS

Model	MC-D06
Overall size	1920*920*600-890mm
Back section up/down	0-65°±5°
Foot section up/down	0-40°
Height adjustment	600-890mm
Trendelenburg	0-10°±2°
Revere- Trendelenburg	0-10°±2°
Foot section fold up/down	55°±2°/5°±2°, removable
Safe working load	200kg
N.W	100kg

ACCESSORIES

Waterproof cushion

Waist grab handle

Leg holder

Filth basin

Clamp

I.V. pole

Power cord

MAIN DETAILS



Leg section adjustment

The middle part of the leg section can be pushed under the bed and quickly adjusted to the position for delivery and gynecological examination



Filth basin

Removable full stainless steel waste basin with drain pipe



Foldable leg section

The leg section can be folded up into footrest and can swing out and removable



Leg rest adjustment

Height of the leg rest can be adjusted, with a universal ball, the angle can be adjusted in all directions



Hand grip

Height of hand grip can be adjusted



Side rail

Two split side rails, tuck-away type



Angle display

Angle display on the side rail, can indicate the adjustment angle of the bed.



Control Panel

Guardrail with integrated control panel



Bumper castors

Bumper castors on the head of the bed



Removable headboard

Headboard is removable and lockable



Castors

Four 125mm three-position central brake castors, with guidance function



CPR function

Manual CPR handle on both the side of the back section. Release the back-rest for emergency situations.



Seat section adjustment

Seat section can be adjusted from 0-12° by remote



Address : No.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng
Town,Zhangjiagang,Jiangsu Province,China
Tel :139 6223 1292
E-mail : info@medik-medical.com