

MC-D08

Gynecological Table



*Better care
Better comfort!*

MC-D08

GYNECOLOGICAL TABLE

DESCRIPTION

MC-D08 is a gynecological table used for delivery and gynecological procedure.

Compact design with electrically adjustable height provides doctor great convenience to get closer operation. It is easy to move the table to designed positions with big lockable casters. It is an ideal option for gynecology option with limited budget.

MAIN FEATURES

- The Structure is epoxy coating
- Mattress adopts one-time foam molding, easy to clean and sterilize.
- Back, seat and height adjustable by foot pedal(optional)
- Stainless steel waste bin
- Four swivel casters with locking casters in rear
- Leg Supports and Side Loops fitted as standard

TECHNICAL PARAMETER

Model	MC-D08
Size	1450*600*600-1000mm
Backrest adjustable	-17° to 70°
Trendelenburg	17°
Seat section adjustable	18°
Load capacity	180KG

PRODUCT DETAILS



Mattress

The mattress adopts one-time foam molding, easy to clean and with beautiful outlook.



Hand grip

The hand grip can be adjusted according to the needs of pregnant women, with ribbed for added resistance.



Foot support

The foot support is easy switch to be positioned to various positions.



Dirt basin

The dirt basin can be stretched and adjusted to different positions, made of stainless steel, easy to clean.



Castor

Four swivel casters make it easy to move to the desired position.



Foot plate (optional)

There is plate on the legs as an option, which can be removed to meet different needs.



Foot control (optional)

Foot control function as an optional choice, to satisfy different doctors' needs.



Roll holder (optional)

Optional roll holder make it more easy to change the bed sheet.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.