

# MK-A06

Hospital Accompanying  
Sleeping Chair



*Better care  
Better comfort!*

# MK-A06

## HOSPITAL ACCOMPANYING SLEEPING CHAIR

### DESCRIPTION

This accompanying chair has a simple design and can be converted into a simple bed. It has wheels for easy transportation and the overall material is easy to clean, making it an indispensable configuration in the ward.

### FEATURES

- Ergonomic design,with armrest and headrest.
- Can be unfolded into a bed and retracted into a chair.
- One-handed operation, convenient and fast.
- Designed with wheels for easy transportation.

### STANDARD SPECIFICATION

Model	MK-A06
Folded size	D670*W590*H910mm
Unfolded size	L1830*W590*H600mm
Castor diameter	Dia 40mm



**Material and structure**

Steel with epoxy painted frame, foam mattress covered with leather, easy to clean.



**Design**

It can be pulled open to become a companion bed, making it convenient to accompany the patient.



**Castor**

There are 6 castors with rubber feet, which are convenient for deployment, retrieval and transportation, while ensuring stability.



**Pull cord design**

The pull cord design allows you to retract the accompanying chair with one hand.



**Humanized design**

With armrests and headrest, convenient for caregivers and improved safety.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.