

The infusion chair, like any other object in the infusion room, provides the best care and comfort for patients. This MK-F03 infusion chair can adjust the back of recliner and is equipped with strong bottom structure to better protect the safety of the patient.



HANDLE ADJUSTMENT

The chair can be quickly converted into recliner through the lever, and adjustable handle can be converted into chair when you get up. It is convenient for patient to rest during the infusion.



RUBBER FEET

processed epoxy paint and is equipped with 4 pcs of rubber feet.



SPONGE AND LEATHER

High-density sponge with thickness of 4cm, cover leather is flame-retardant, waterproof and wear-resistant.



HAND REST

The wooden hand rest with ergonomic design greatly improves the patient's comfort.

▪ **Size**
68*79*110cm

▪ **Structure**
Epoxy coated steel

▪ **Back section adjustable**
0-90°



Zhangjiagang Medi Medical Equipment Co., Ltd
Tel: +86 15250358182
E-mail: info@medik-medical.com

Infusion chair

Model: YA-MK-F03

- The frame is made of powder painted steel
- With environmental protection PVC leather
- A puller on the right side of chair to control the backrest and legrest simultaneously (gas spring operated) to 3 different positions
- With a stainless steel I.V pole with 2 hooks

