

MK-F05

Infusion chair



*Better care
Better comfort!*

MK-F05

Infusion chair

DESCRIPTION

MK-F05 is a luxurious blood donor chair as well as for infusion chair with swing away arms , the special design provides both patients and nurses with comfortable experience specially for disabled or extra size patient. Designed specifically for the clinical environments, including uses in Dialysis, Oncology & Infusion, Orthopedics, Same-day Surgery, Long-term Care, IV Therapy, PACU and Emergency Room.

FEATURES

- The puller on the right side of the chair can swing away arms
- A puller on the right side of chair to control the backrest and legrest simultaneously to 3 different positions
- Push handle to move the chair
- With a stainless steel I.V rope with 4 or 2 hooks

STANDARD SPECIFICATION

Overall size	1280*850*880mm
Seat size	500*570*540mm
Bed size	1750*570mm
Backrest adjustable	0-75°
Footrest adjustable	-75- 0°
Load capacity	150 KG



BACK AND FOOT SECTION ADJUST

The puller on the right side control the angle of back and foot section , Toggle the handle to control the back and legs and adjust it to a comfortable angle for the user.It can also adjust as a bed position.



SWING AWAY ARMS

The puller on the left side can swing away the arms , Toggle the handle to move away the left arms which is more convient for patient to get on and off the chair .



SIDE DINING TABLE

Dual side dining table self-store by folding down and are easily cleaned.



UPHOLSTERED ARMRESTS

Upholstered Armrests in matching fabric offer clean lines and attractive.Have a hole to set the bottle or something else in case it is accidentally knocked off.



PUSH HAND

The push hand behind the back section can better enable the staff to move the chair.



IV STAND

The height-adjustable IV stand can meet the needs of different heights patients.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US