

YA-DS-D02

Medical Hospital Dialysis Chair



*Better care
Better comfort!*

YA-DS-D02

Medical Hospital Dialysis Chair



Our electric blood donor chair YA-DS-D02 is designed for blood donation/collection, blood transfusion or intravenous (IV) treatment.

TECHNICAL PARAMETER

Model	YA-DS-D02
Overall size	1980*550*580-750mm
Back section adjustment	0~ 70°
Leg section adjustment	20° -90°
Height adjustment	580-750mm
Armrest swing out	0 °-180°
Armrest height adjustment	100mm
Armrest size	480*160mm
Loading capacity	180kg

- Ergonomic design gives blood donors a sense of security and comfort and alleviate the psychological pressure, makes the whole process easier and enjoyable.
- The appearance is made of high-strength fiber material and poly-layer paint
- Upholstery with 60mm thick foam
- Equipped with 3 motors for height, back and leg rest adjustment
- Chair frame is made of premium powder coated steel
- Armrests adjustable in height, depth and with pivoting articulation
- Four \varnothing 100 mm casters with individual locking system are included to facilitate the movement of the donor chair.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US