

YA-DS-M02

Infusion chair



*Better care
Better comfort!*

INFUSION CHAIR

YA-DS-M02

YA-DS-M02 is a luxurious infusion sofa with one hand back adjusting. The ergonomic sofa lines can bring the user a comfortable treatment experience to the greatest extent. Caster is optional for easy reposition



Main Features

- Etherealized lines, freely outline the casual, fashionable and graceful form
- Environment-friendly and breathable material, anti-fouling and dust proof anti-static, new for long lasting
- Material: Solid wooden frame with multi-layer iron plate
- High resilient cutting foam and PU cover, thick foam bring comfort
- One-handed easy adjustment of body position, convenient and comfortable
- With IV pole

Technical parameters

Model	YA-DS-M02
Overall size	1750*820*1050mm
Folding size	950*820*1050mm
Back section adjustment	0°-80°
Caster	Diameter 50mm
Loading capacity	200kg

Main details



Leg and back section is adjusted by mechanism adjustment switch



Mechanism adjustment switch




Cup holder



Thick foam upholstery for comfort



Bottom can be castor or plastic feet



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US