

YA-DS-M09

Dialysis Recliner Chair



*Better care
Better comfort!*

YA-DS-M09

DIALYSIS RECLINER CHAIR

DESCRIPTION

YA-DS-M09 is an dialysis recliner chair with ergonomic design for improving patient comfort. Plus the high density cushioned mattress and wider armrest with cup holder and practical foldable table, it is a very popular chair for the dairy dialysis (or infusion) treatment in the clinic.

FEATURES

- Epoxy coated durable steel frame structure.
- Upholstered with High-Quality Leather & Padded with 8cm High-Density Foam - cover water proof, easy cleaning and sterile.
- For easy Move & Transport provided with High Quality 4 castors and brake.
- Self control adjusted up the feet & back resting with a gas pressure spring mechanism place to provide a safe & comfortable environment for blood donation or for Patient Assistant .
- Seat Size: W58cm x D54cm.
- Unfolding Size: 170x80x45cm.

STANDARD SPECIFICATION

Overall size	90*85*105cm
Seat size	58*54cm
Back-leg adjustment	90-165°
Loading capacity	180kg



Handle

The adjustment of the back and legs can be adjusted by manual handle.



Small dining table

The dining table is on the right side of the chair, it can place small items and be folded.



Water cup

The cup is embedded in the left side of the armrest of the chair, where you can place a water cup.




IV pole holder

The IV pole holder is behind the chair, used to place the infusion pole.



Castor

Dia 100mm castors , 2 with brakes.



IV pole

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.