

HOSPITAL NURSE STOOL

Product introduction

YA-S01 is a star base structured stool with foot rest, excellent mobility with five casters. The height is conveniently adjusted in hydraulic way. It is an ideal moving stool for doctor checking room, nursing station or lab.

Technical parameter

Size	Ø35*45-61cm
Material	Chrome plated structure with PU foam seat top
Castors	Five 50mm castors
Height adjustable	45-61cm by gas spring



info@medikmedical.com | +86 13962231292

HOSPITAL NURSE STOOL

Main details



Main material

The main material is the high quality stainless steel, which is easier to clean, modern and stable.



Seat surface

PU foam seat top, round, beautiful and can meet any need, easy to clean.



Height adjustment

Height adjustment can be from 45-61cm to meet different height needs. Safe and convenient control by push handle.



Basement

The bottom is rounded, which can rest the feet, improve comfort and high stability.



Castor

Dia swivel 50mm castors which move flexibly.



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table&lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.



YOU TUBE MEDIK

Pls scan the QR code for more info

info@medikmedical.com | +86 13962231292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG - US