

YA-PS17

Foldable Patient Transfer Stretcher



YA-PS17 FOLDABLE PATIENT TRANSFER STRETCHER

The YA-PS17 Foldable Patient Transfer Stretcher is designed with a space-saving folding structure, allowing hospitals and healthcare facilities to easily store and deploy the trolley when needed. Featuring a sturdy integrated frame, adjustable backrest, rotating armrests and smooth-rolling casters, it provides safe and comfortable patient transfer while improving storage efficiency in limited spaces. With a built-in IV stand holder and reliable braking system, this foldable transfer trolley is ideal for emergency rooms, clinics, wards and other medical environments where flexible patient mobility is required.

Main features

- Foldable space-saving frame design for easy storage and transport
- Heavy-duty integrated frame for reliable patient transfer
- Rotatable height adjustable lockable armrests
- Manual backrest adjustment for comfortable positioning
- Integrated IV pole holder for infusion support
- Four silent casters with individual locking brakes

Main details



Frame

Collapsible integrated metal frame, easy unfolding & folding for space-saving storage.



Armrest

Rotatable height-adjustable armrests with bayonet locking for firm fixation.



Backrest

Manually smooth adjustable backrest to meet different lying postures.



Castors

Four silent casters with independent brakes, stable and anti-slip during use.



Accessory

Reserved IV stand holder, compatible with infusion bracket installation.

Technical parameter

Model	YA-PS17
Unfolded Size(LWH)	1880×700×640mm (74" x 27.6" x 25.2")
Folded Size(LWH)	680×350×1140mm (26.8" x 13.8" x 44.9")
Castors	Dia120 mm(5")
Backrest Adjustment	0-18°



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US