

MK-S43

Stainless Steel Instrument Trolley



A LEADING BRAND FOR HOSPITAL FURNITURE SOLUTION

+86 139 6223 1292 | www.medik-medical.com

medik⁺

MK-S21 STAINLESS STEEL INSTRUMENT TROLLEY

The MK-S43 Stainless Steel Instrument Trolley is designed for hospitals, clinics and operating rooms, providing a practical solution for transporting and organizing medical instruments. Constructed from premium 304 stainless steel, it offers excellent corrosion resistance, easy cleaning and long-lasting durability for daily clinical use.

Main features

- Premium 304 stainless steel construction.
- Two stainless steel shelves for efficient instrument storage.
- Raised tray edges with three-sided guard rails help prevent items from slipping.
- Ergonomic integrated push handle for easy transportation.
- Smooth brushed surface with polished welds for easy cleaning.
- Four 100mm medical silent castors, two with brakes.
- Stable frame construction suitable for daily hospital use.
- Customization options available upon request.

Main details



Premium Stainless Steel Construction

Made of high-quality 304 stainless steel with a brushed finish, providing excellent corrosion resistance, easy maintenance and a professional appearance.



Double Shelf Design

Two spacious stainless steel shelves provide sufficient storage space for medical instruments and supplies, improving working efficiency.



Raised Tray with Guard Rails

Both shelves feature raised edges and three-sided guard rails to help keep instruments securely in place during transportation.



Integrated Push Handle

The ergonomic push handle allows smooth and comfortable maneuvering while maintaining excellent trolley stability.



Medical Silent Castors

Equipped with four 100mm medical-grade swivel castors, including two locking castors, ensuring quiet movement and secure positioning.

Technical parameter

Model	MK-S43
Material	Stainless Steel 304
Overall Size	600×400×800 mm (23.6" × 15.7" × 31.5")
Shelves	2 pcs
Castors	Dia100 mm (4")
Brakes	2 pcs
Handle	Integrated Push Handle
Application	Hospital / Clinic / Operating Room



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US