

MK-BH01

Medical Infant Height Bed



MK-BH01 MEDICAL INFANT HEIGHT BED

The MK-BH01 features a robust spray plastic steel construction, which is smooth, rust-proof and durable for long-term service. It provides a 0-100cm measuring range with 0.1cm high-precision scale, adopting a stable and easy-to-use design ideal for height measurement of 0-3-year-old infants in hospitals and clinics.

Main features


- Sturdy spray plastic steel frame, rust-proof and durable
- Smooth edge design, safe for infant use
- Clear and high-precision measuring scale
- 0-100 cm measuring range with 0.1 cm accuracy
- table structure for comfortable lying
- Simple and labor-saving operation

Main details



Frame

Spray plastic steel frame, rust-proof, sturdy and deformation-resistant.



Surface

Smooth burr-free surface, safe for infant skin contact.



Division

High-precision clear division, minimum value of height per division by 0.1cm.



Measuring Range

0-100cm measuring range for 0-3 years old infants.

Technical parameter

Model	MK-BH01
Measuring Range	0-100 cm
Material	Spray Plastic Steel
Minimum Scale Value	0.1 cm
Applicable Age	0-3 Years Old



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US