

YA-S04

Multi-Function Doctor Chair



A LEADING BRAND FOR HOSPITAL FURNITURE SOLUTION

+86 139 6223 1292 | www.medik-medical.com

medik⁺

YA-S04 MULTI-FUNCTION DOCTOR CHAIR

The YA-S04 ergonomic multi-function doctor chair serves as auxiliary passive equipment for clinical examination and diagnosis. It features cross-shaped backrest and movable saddle-shaped seat to ease operators' physical fatigue. With hydraulic gas lift, multi-control base and locking casters, it enables multi-directional adjustment, helping medical staff conduct operations from various angles. Suitable for examination rooms and clinical diagnostic environments.

Main features

- Human-centric ergonomic design, cross-shaped backrest and movable saddle-shaped seat
- Seat cushion fore-aft sliding: 0-60 mm; backrest tilting: 0-10°, with lock & auto-reset function
- Hydraulic seat height: 510-700 mm; adjustable backrest height: 920-1120 mm
- Multi-functional manual-control base for convenient position adjustment
- Steel frame treated with double electrostatic anti-rust and antibacterial coating
- Heavy-duty five-star aluminum-alloy base (365 mm radius), superior stability
- Certified durable PU upholstery, high-density sponge inner filling
- Heavy-duty gas lift rod, passed 10000-cycle service life test
- 2-inch silent casters with seat-brake locking function
- Maximum safe load capacity: ≥170 kg



Main details

Ergonomic Seat & Backrest

Adopts cross-shaped backrest and saddle-shaped seat cushion. The seat slides forward and backward within 0-60 mm; backrest tilts from 0-10° with position locking and automatic reset upon unlocking. It effectively relieves lumbar strain and shoulder discomfort for long-hour working medical personnel.



Height Adjustment

Hydraulic seat height ranges 510-700 mm. Backrest height adjusts from 920 mm to 1190 mm, operated via multi-functional manual-control base for quick posture change.



Heavy Frame & Five-Star Base

High-quality steel frame with double electrostatic coating for anti-rust, antibacterial and anti-corrosion performance. Heavy-duty 5-spoke aluminum-alloy base (radius 365 mm) delivers solid, stable support.



Upholstery

High-density shaped sponge padding, covered with certified water-based solvent-free PU synthetic leather for comfort and wear resistance.



Hydraulic Gas Lift Rod

Reliable gas lift rod with ≥ 100 Kg load capacity, passed 10000-times lifting cycle durability test.



Silent Brake Casters

2-inch deluxe silent casters equipped with seat-activated lock, preventing undesired sliding during diagnosis for higher safety.



Technical parameter

Model	Multi-Function Doctor Chair
Overall dimension	730x730x920-1120 mm
Seat cushion fore-aft travel	0-60 mm
Seat height adjustment	510-700 mm
Backrest tilting angle	0-10°
Backrest overall height	920-1190 mm
Max width of side wings on backrest	540 mm
Seat dimension	Width 460 mm, Depth 400 mm
Base	Heavy-duty Five-star aluminum-alloy base, radius 365 mm
Caster	2-inch deluxe silent caster with seat-brake lock
Gas rod load	≥100 Kg, 10000 lifting cycles tested
Overall safe load	≥170 Kg
Padding	High-density shaped sponge
Surface upholstery	Certified water-based solvent-free PU synthetic leather
Frame finish	Double electrostatic coating, anti-rust, antibacterial, anti-mildew, anti-corrosion



Unit: mm

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 30 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia. Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval. Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US