

YA-EC-U01

Electric Ultrasound Examination Table



*Better care
Better comfort!*

YA-EC-U01

ELECTRIC ULTRASOUND EXAMINATION TABLE

DESCRIPTION

Automatic linen sheet change will improve exam efficiency, especially for the ultrasonic room, and it will provide a good hygiene experience for patients. YA-EC-U01 is an ideal electric examination couch for ultrasonic room with big storage space and an automatic linen change system.

MAIN FEATURES

- Steel with powder coating structure, durable, antibacterial and anti-corrosion
- Seamless high-density mattress, antibacterial,corrosion-resistant and scratch-proof
- With huge storage space
- Sheet roll holder on both ends of the bed
- Automatic bed sheet change with automatic deviation correction function
- Quick bed sheet changing system
- Four castors ,dia 5 inches ,for easy transport
- With battery backup for emergency situation

TECHNICAL PARAMETER

Model	YA-EC-U01
Overall size	2010*720*620mm
Static working load	250kg
Dynamic working load	170kg
Castor	Dia 5 inchies
N.W	83kg

PRODUCT DETAILS

Pillow and Mattress

Ergonomically designed pillow for easy use,seamless high-density mattress,antibacterial,corrosion-resistant and scratch-proof ,easy to clean.




Push handle

Push handle on the end of the bed ,easy to transport.



Sheet roll holder

Sheet roll holder on both ends of the bed ,easy for automatically sheet change.



Quick bed sheet changing device

There are quick bed sheet changing devices , making it easy to replace new sheets.



PRODUCT DETAILS

Hand remote

The hand control is used to automatically change sheets. It has two speed, has an automatic correction function, and is equipped with a disposable mattress.



Storage space

Lockable storage space, can be opened from both sides.



Castors

Four pcs of dia 5 inches castors with individually brake.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US