

medik



MC-C01

GYNECOLOGY
EXAM TABLE

BETTER CARE
BETTER COMFORT

Zhangjiagang
Medi medical Equipment
Co., Ltd



You Tube



Medik

Tel : +86 15250358182 / Fax : +86 51258688450

E-mail : info@medik-medical.com

Address : No.3 Jinxiu Road,Jinfeng Industrial Zone,Jinfeng

Town,Zhangjiagang,Jiangsu Province,China

MC-C01 / Gynecology exam table

The MC-C01 gynecological exam table manages the room space maximally by integrated cabinet and extensible leg platform and foot support, it adopts a seamless design to reduce the accumulation of dirt and bacteria. This Patient Exam Table maximizes convenient features including built-in stirrups, footstep that can be easily pulled out or tucked away. This table is a perfect option for doctor's offices, clinics or wherever a cost effective treatment or exam table is needed. .



SIDE DRAWERS

Three side drawers provide point of care storage for gowns, drapes and other items. It has 4 casters and can move as a whole.



DIRT BASIN

The pull-out dirt basin is very suitable for lower body operation that need to remove liquid and debris.



FOOTSTEP

The bottom of the front of MC-C01 is equipped with pull-out footstep, and the surface of the footstep is non-slip treatment, which can make the whole more safe and convenient.



BACK SECTION ADJUSTMENT

The back section is adjusted by gas spring, which can be adjusted from 0-70°.



STIRRUPS

Pull out stirrups provide leg support for patients during lower body examinations.



PAPER ROLL HOLDER

The paper roll holder is hidden behind the back section and can be seen by lifting back section.

NO.	ITEM	PARAMETERS
01	Overall size	1178*520*520mm
02	Back section adjustment	0-70°
03	Seat section adjustment	0-7°
04	Step size	420*230mm(L*H)
05	Load capacity	180kg

MAIN FEATURES

- The ob gyn exam table is suitable for all kinds of inspection including gynecological examination, the shape is beautiful, the function is practical, the operation is simple.
- The back section is operated by gas spring.
- The roll paper holder is hidden behind the back panel.
- Foot stool can be used for patient as a step, which is hidden.
- Examination bed is equipped with large drawers.
- Mattress made from mould, with no seam on the surface.