

# MC-C02

## Gynecological Examination Table



*Better care  
Better comfort!*

# MC-C02

## Gynecological Examination Table



### DESCRIPTION


The MC-C02 examination table has a large storage capacity and with necessary accessories for gynecological procedure. With retractable auxiliary leg section and foot step , It is a perfect option for doctor’s offices, clinics or wherever a cost-effective treatment or exam table is needed.

### STANDARD SPECIFICATION

Overall Size	1870*600*680mm
Packing Dimensions	1520*710*830mm
Back-rest lifting angle	75° (±5°)
Gross Weight	92 kg
Loading capacity	150kgs

**BACK SECTION ADJUSTMENT**

Back section can be adjusted by ratchet system from 0° to 70°.



**FOOT SECTION**

Foot section can be hidden under the seat section, unfolded to raise the foot section as high as the seat section.



**STRUCTURE**

Base cabinet structure, Electrophoresis and powder coating, thickness 1 mm to ensure the stability of the work table. load capacity 250kg.



**PILLOW**

The pillow has a unique shape and is more suitable for the patient. the pillow is fixed with bristles and can be taken freely.



**FOOT STEP**

The upper and lower exam beds of the patient can be accessed by the foot step, and it also with anti-slip pad.



**DRAWERS**

The gynecological examination table has large capacity drawers that can be used to place inspection instruments and documents.



**LEG SUPPORT**

The gynecological examination table is equipped with 2pcs leg support, which can be used for some routine examinations.



**MATTRESS**

Modular foaming high-density mattress makes people more comfortable, easy to hygiene.



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US