



Tel : +86 15250358182
E-mail : info@medik-medical.com

SIMPLE GYNECOLOGICAL EXAMINATION CHAIR MK-C06

Simple gynecological examination chair, using high-quality stainless steel plate and round tube with argon arc welding, smooth and bright, without concave deformation, can meet the needs of most gynecological examinations.

Leg section adjustment

The leg section can be adjusted vertically 90 degree down so the doctor can get closer with convenience, it also can be adjusted in horizontal position with seat section.



Leg support



The leg support can be adjusted in height and curvature, which greatly improves the comfort of the patient.

Dirt Basin



Stainless steel dirt basin easy to clean, durable and replaceable.

Mattress



3cm thick mattress, the surface is wrapped in PVC leather, easy to clean.

Handle



The curved handle can adjust the height and provide support during the patient's delivery or examination.

Bed legs



Rubber legs are installed on the bed legs to prevent noise from rubbing with the floor.

- Overall size : 1900*610*750mm
- Back section adjustment : 0-75°
- Leg section adjustment : 0-90°
- Loading capacity : 170kg