

MC-D11

Gynecologist Chair



*Better care
Better comfort!*

MC-D11

GYNECOLOGIST CHAIR

DESCRIPTION

The MC-D11 is a gynecologist chair designed for comprehensive health examinations. With versatile features and several control ways, it promises not only beauty but also unparalleled performance. this chair ensures a comfortable and efficient experience for patients and healthcare professionals.

FEATURES

1. Featuring a built-in foot pedal, the MC-D11 frees hands, minimizing hygiene concerns.
2. The one-key reset function facilitates a swift return to the initial state, allowing patients to quickly exit the chair
3. Enhancing safety measures, the lock button prevents accidents resulting from inadvertent foot pedal activation, ensuring a secure environment.
4. Armrests linkage with broader applicability, catering to various sitting, lying down, or front-lying positions.

STANDARD SPECIFICATION

Size	1670*825*570-895mm
Backrest adjustable	0 to 80°
Footrest adjustable	0 to 90°
Seat section adjustable	25°
Load capacity	200KG



Built-in foot pedal

Featuring a built-in foot pedal which empowers operators by freeing their hands, thereby minimizing hygiene concerns.



Headrest adjustable

The headrest pillow is adjustable and designed to cater to a diverse range of needs and preferences.



Arm rest

The armrest flips upward, facilitating easy patient exit from the table.



Dirt basin

The dirt basin can be stretched and adjusted to different positions, made of stainless steel, easy to clean.





Foot support

The foot support is easy switch to be positioned to various positions .



Leg adjustable

The leg section adjustment is controlled by gas springs, ensuring simple, convenient, and swift operation.



Back Cover

The meticulously crafted backrest cover not only enhances the bed's visual appeal but also offers practicality in terms of easy cleaning and disinfection.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.