

MC-D18

Gynecological Chair



*Better care
Better comfort!*

MC-D18

GYNECOLOGICAL CHAIR

DESCRIPTION

The MC-D18 is a gynecological chair designed for birthing delivery, gynecological exam. The unique recessed frame structure and foot switch control provides the doctor a great convenience to do the operation.

FEATURES

- The Structure is epoxy coating
- Seat section folds up by gas spring.
- Back and height adjustable by foot pedal
- Stainless steel 5L waste bin
- Four swivel casters two with brakes
- Leg Supports and Side Loops fitted as standard

STANDARD SPECIFICATION

Size	1460*750*550mm
Height adjustable	550-1050mm
Backrest adjustable	-20° to 70°
Trendelenburg	20°
Seat section adjustable	30°
Castors	100mm
Load capacity	180KG



Arm rest

The armrest is meticulously designed to not only provide support but also enhance overall comfort, ensuring a relaxed and comfortable experience for patients during examinations.



Leg support

The leg support is easy switch to be positioned to various positions.



Dirt basin

The dirt basin can be stretched and adjusted to different positions, made of stainless steel, easy to clean.



Seat plate

The seat plate is adjusted by the gas spring and seamlessly adjustable at various angles, providing flexibility for tailored patient positioning during examinations.



Foot rest

Footrests are meticulously designed to deliver enhanced patient comfort, ensuring optimal relaxation during examinations.



Foot control

With the foot-operated height adjustment which is intuitive control, allows doctors to focus on the procedure of checking the patient care.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US